

NDC 43386-357-10

**Hydrocodone Bitartrate
and Acetaminophen
Tablets, USP**



7.5 mg/325 mg

Pharmacist: Please dispense
accompanying Medication Guide
to each patient.

Rx only
LUPIN

1000 Tablets

Each tablet contains:
Hydrocodone bitartrate, USP.....7.5 mg
Acetaminophen, USP..... 325 mg

Pharmacist: Dispense in a tight, light-resistant container with a child-resistant closure.

Usual Adult Dosage: One tablet every four to six hours, as needed for pain. Total daily dosage should not exceed six tablets.

See insert for full prescribing information.

Keep this and all medications out of the reach of children.

Storage: Store at 20° to 25°C (68° to 77°F) [See USP Controlled Room Temperature].

Manufactured by:
Novel Laboratories, Inc.
Somerset, NJ 08873

Manufactured for:
Lupin Pharmaceuticals, Inc.
Baltimore, MD 21202

LA3571000201

Iss. 06/2017



38mm x 153mm
Unvarnished Area



CONTROL GROUP®

Printed Packaging Solutions: Pharmaceutical • Cosmetic • Consumer
500 Walnut Street, Norwood, N.J. 07648 • www.controlgroupusa.com

proof date	06/14/2017
proof time	10:37 AM
operator	maureen.pedersen

Customer: **Lupin - Hydrocodone Bitart & Acet Tabs 7.5mg-325mg X 1000**

P.O. No #: **122185** Cust. Part#: **LA3571000201**

Job number: **NOVE_0304** Work Order: **35950**


Size: **4.000 x 11.000 in (101.60 x 279.40 mm)**

ELECTRONICALLY SPELL CHECKED BY:	PROOFREAD INTERNALLY BY:	DATE:
-------------------------------------	--------------------------	-------

Reviewed by:	Date
Submit Revised Proof:	Date
Proof Approved By:	Date

Inks:


NOTE: Proof colors do not represent exact PMS colors. Please check current PMS guide.



Non-printing technical colors: **Redline** **Unv Area**

Web:

Roll direction chart



Roll Direction #4

Approved for:

	Yes	No	N/A
1 Copy	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
2 Web layout (inc. eyemark)	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
3 Roll directions	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>